



Zhejiang AnyGo Technology Co., Ltd.
107-1, Area A, Building 1, No.6 Kaixuan South Road, Green
Industry Cluster District, Quzhou City, Zhejiang, China.

AK200-Mouse Rabies Virus Protective Antibody Detection Kit Instructions for AK200 ELISA Kit

【Purpose】

The AK200 ELISA Kit is designed for semi-quantitative detection of rabies virus protective antibody in mouse serum samples, for in vitro research use only.

Cat# AK200M1

【Experimental principle】

Briefly, flat-bottom 96-well plates are coated with Glycoprotein of rabies virus that captures rabies virus protective antibody in the sample. A quick wash removes any unbound particles. Captured rabies virus protective antibody is detected by Goat Anti-Mouse IgG pAb, which is conjugated to horseradish peroxidase (HRP). Finally, the chromogenic substrate 3, 3', 5, 5'-tetramethylbenzidine (TMB) is added. The amount of rabies virus protective antibody is proportional to the color generated in the coupled oxidation-reduction reaction compared with known amounts of standards.

【Kit Components】

96-well plate coated with Glycoprotein, 8 wells×12 strips

Negative control, 40μL/tube

Positive control, 40μL/tube

Sample Dilution Buffer, 12mL

Goat-anti-Mouse IgG pAb (HRP Conjugated pAb, 100X), 120μL

BSA, 3g

20X PBS Buffer, 50mL

Substrate Solution A, 6mL

Substrate Solution B, 6mL

Stop buffer, 6mL

Sealing films, 2 pieces

Assay Instruction

【Storage Conditions】

When stored at 2~8°C, the product is stable for one year.

【Protocol】

- A. Before using, pre-warm all the reagents to room temperature (18~25°C).
- B. **Washing buffer** dilution: dilute 20X PBS Buffer to 1000mL and mix well.
- C. **Enzyme Dilution buffer (for HRP conjugated pAb)**: dissolve 3g BSA to 100mL washing buffer and mix well.
- D. Dilute samples with **Sample Dilution Buffer** and add 100µL into each well of the ELISA plate. Set negative control, positive control and incubate at 37°C, 200 rpm for 60min.
- E. Remove samples from wells and wash all wells three times with prepared Washing Buffer. Remove residual solution with paper pat.
- F. Dilute 100µL Goat-anti-Mouse IgG pAb with 10mL **Enzyme Dilution buffer (for HRP conjugated pAb)** and mix well. Add 100µL of diluted Goat-anti-Mouse IgG pAb to each well. (Careful not to touch or scratch the surface of the wells). Incubate plate at 37°C, 200rpm for 60min.
- G. Remove samples from wells and wash all wells three times with prepared Washing Buffer. Remove residual solution with paper pat.
- H. Add 50µL each of Substrate solution A and Substrate solution B into each well. Mix thoroughly with shaking. Incubate at 25°C for 5 min and avoid light.
- I. Stop the reaction by adding 50µL of Stop Solution.
- J. Record the absorbance at 450nm on a plate reader within 30min.

【Properties】

1. Precision: coefficient of variation CV (Percentage) should not be higher than 15 % (n=10).
2. Stability: The components of the kit are left for 7 days at 37°C under accelerated destruction condition, and all performance indexes still meet the above requirements.

【Caution】

- A. Avoid cross contamination.
- B. Follow reader measure as standard.
- C. All samples and buffers are added or removed with pipette.
- D. Do not mix reagents from different batches.