

## Mouse anti-Meningococcal polysaccharide(Y) antibody, mAb 7C6

**Cat#: 05-0103**

GENERAL INFORMATION	
Physical State	Liquid
Immunogen	Meningococcal polysaccharide(Y)
Preparation	This antibody is produced by a hybridoma cell, which is produced by the fusion of a mouse myeloma and a B cell, the B cell is obtained by immunizing mice with Meningococcal polysaccharide(Y). The IgG fraction of the cell culture supernatant was purified by protein G affinity chromatography.
Specificity	Base on ELISA, this antibody reacts with Meningococcal polysaccharide(Y) Specifically
Formulation	Please see Certificate of Analysis for specific instructions. For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.
Storage	This antibody can be stored at 2°C-8°C for short term (1-2 Week) without detectable loss of activity. Antibody products are stable for 2 years from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications	ELISA
RECOMMENDED CONCENTRATION	
ELISA	This antibody is suitable for use in ELISA: $\geq 0.5\mu\text{g/mL}$ .
Please Note: Optimal concentrations/dilutions should be determined by the end user.	

AbMax Biotechnology Co., Ltd 99 Kechuang 14<sup>th</sup> St, Bldg 18, 2-201 BDA, Beijing, China 101111

[www.antibodychina.com](http://www.antibodychina.com)

Not for use in diagnostic procedures. If you have any further questions, please contact Technical Services at 010-59755726-8006

**For Research Use Only.**(仅供体外科研使用)