

Anti-Varicella zoster virus(VZV) mAb 5E6,HRP-Conjugate

Cat#: 12-0115

GENERAL INFORMATION	
Physical State	Liquid
Description	HRP Conjugate Mouse anti-Varicella zoster virus(VZV) antibody, mAb 5E6.
Specificity	Base on ELISA, this antibody reacts with gE Protein Specially.
Conjugate	HRP
Formulation	Please see Certificate of Analysis for specific instructions(0.01M PBS. pH7.2, 0.01% sodium azide, 20%Glycerol). For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.
Storage	This antibody can be stored at 2°C-8°C for short term (1-2week) without detectable loss of activity. Antibody products are stable for 2 years from date of receipt when stored at ≤ -15°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications	ELISA
Please Note: Optimal concentrations/dilutions should be determined by the end user.	

AbMax Biotechnology Co., Ltd 99 Kechuang 14th St, Bldg 18, 2-201 BDA, Beijing, China 101111
www.antibodychina.com

Not for use in diagnostic procedures. If you have any further questions,
please contact Technical Services at 010-59755726-8006

For Research Use Only.(仅供体外科研使用)