

Anti-varicella zoster virus(VZV) antibody, mAb 5E6 ,Biotin-Conjugated

Cat#: 12-0109

GENERAL INFORMATION	
Physical State	Liquid
Description	Biotin Conjugate Mouse anti Varicella zoster virus(VZV) antibody, mAb 5E6
Specificity	Base on ELISA, this antibody reacts with gE-FC protein
Conjugate	Biotin
Formulation	Please see Certificate of Analysis for specific instructions(0.01M PBS,pH 7.2,0.01% sodium azide, 20% glycerol). For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.
Storage	This antibody can be stored at 2°C-8°C for short term (1-2week) without detectable loss of activity. Antibody products are stable for 2 years from date of receipt when stored at ≤ -15°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications	ELISA
Please Note: Optimal concentrations/dilutions should be determined by the end user.	

AbMax Biotechnology Co., Ltd 99 Kechuang 14th St, Bldg 18, 2-201 BDA, Beijing, China 101111

www.antibodychina.com

Not for use in diagnostic procedures. If you have any further questions, please contact Technical Services at 010-59755726-8006

For Research Use Only.(仅供体外科研使用)