

Mouse anti-Meningococcal polysaccharide(A) monoclonal antibody

Non WHO Reference Material

Instructions for use (Version 1.0, Dated 2025/3/6)

Catalog No.: 21-0016

This material is not for in vitro diagnostic use

INTENDED USE

Mouse Monoclonal Antibody for use in the Immunoturbidimetric assay for measuring the concentration of Meningococcal polysaccharide(A) .

CAUTION

This preparation is not for administration to humans or animals in the human food chain. This reference contains material from murine origin. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.

IMMUNOGEN

Meningococcal polysaccharide(A)

MATERIALS

Mouse IgG, produced by a monoclonal hybridoma cell and purified from mouse ascites by Protein G column, in 10 mM Tris-glycine-HCL, 0.01% NaN₃.

UNITAGE

1mL/vial.

USE OF MATERIAL

The material can be diluted at 1:2~1:4. However where possible the user should determine the working dilution depending on the requirements of the Immunoturbidimetric.

Example Data:

The material was diluted 1:2, then take 20 ul of the diluted antibody solution to mix with 20 ul of sample and 160 ul of reaction buffer. Gain: 3, Reaction time:2.5 min.

CONTENTS

Country of origin of biological material: People's Republic of China.

STORAGE

The material should be stored frozen at $\leq -15^{\circ}\text{C}$. Once thawed, keep at $2-8^{\circ}\text{C}$ to avoid freezing and thawing. Material type: Liquid – will be shipped according to the storage and shipping conditions of the product

STABILITY

The materials are stable for 2 years at $\leq -15^{\circ}\text{C}$, 1 month at $2-8^{\circ}\text{C}$ and 2 weeks at room temperature.

FURTHER INFORMATION

info@antibodychina.com

TECHNICAL SUPPORT

jszc@antibodychina.com

MATERIAL SAFETY SHEET

| Physical and Chemical properties | |
|---|---|
| Physical appearance: Frozen liquid | Corrosive: No |
| Stable: Yes | Oxidising: No |
| Hygroscopic: Yes | Irritant: No |
| Flammable: No | Handling: See Section 2 |
| Others: Containing 0.01% NaN_3 | |
| Toxicological properties | |
| Effects of inhalation: | Not established, avoid inhalation |
| Effects of ingestion: | Not established, avoid ingestion |
| Effects of skin absorption: | Not established, avoid contact with skin |
| Suggested First Aid | |
| Inhalation: | Seek medical advice |
| Ingestion: | Seek medical advice |
| Contact with eyes: | Wash with copious amounts of water, seek medical advice |
| Contact with skin: | Wash thoroughly with water. |
| Action on Spillage and Method of Disposal | |
| Spillage of ampoule contents should be taken up with absorbent material wetted with an appropriate disinfectant. Rinse are with an appropriate disinfectant followed by water. Absorbent material used to treat spillage should be treated as biological waste. | |

INFORMATION FOR CUSTOMS USE ONLY

| |
|--|
| Country of origin for customs purpose: People's Republic of China, where the goods have been produced and/or sufficiently processed to be classed as originating from the country of supply, for example of a change of state such as freeze-drying. |
|--|

| |
|------------------|
| Net weight: 0.5g |
|------------------|

| |
|---|
| Toxicity Statement: Toxicity not assessed |
|---|

| |
|---|
| Veterinary certificate or other statement if applicable. Attached: No |
|---|