



Anti-MMAE (Monomethyl Auristain E)
(mouse monoclonal IgG)

Catalog # 05-0096B

Immunogen: MMAE

Specificity: Specially recognizes the Monomethyl Auristain E (MMAE)

Formulation: 1mg of mouse monoclonal IgG, purified by Protein G column, in 1ml of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02%NaN₃, 20% glycerol.

Concentration: 1 mg/ml. Protein was determined by OD280

absorbance.

Physical State: Liquid

京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京经济技术开发区科创十四街 99 号 18 号
楼 2 单元 201

99 Kechuang 14th St, Bldg 18, 2-201

BDA, Beijing, China 101111

Tel: 86-10-59755729, Fax: 86-10-59755718

E-mail: info@antibodychina.com

www.abmaxcn.com

Storage: Store the product at -20 °C . Product is stable for about 6 weeks at 2-8 °C as an undiluted liquid. Prepare working dilution prior to use. Avoid repeated freezing and thawing.

FOR IN VITRO RESEARCH USE ONLY

-
- not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user
-