



京天成生物技术（北京）有限公司
AbMax Biotechnology Co., LTD
北京经济技术开发区科创十四街 99 号 18
号楼 2 单元 201
99 Kechuang 14th St, Bldg 18, 2-201
BDA, Beijing, China 101111
Tel: 86-10-59755729, Fax: 86-10-59755718
E-mail: info@antibodychina.com
www.antibodychina.com

Anti-DM1 mAb (1A6)

(mouse monoclonal IgG)

Catalog # 05-0097

Lot #200090

Storage: Store the product at -20 °C .
Product is stable for about 6 weeks at
2-8 °C as an undiluted liquid. Prepare
working dilution fresh each day. Avoid
repeated freezing and thawing.

Formulation: 0.1mg of mouse monoclonal IgG in 0.1ml
of 0.1M Tris-HCl, pH8.0, 0.15M NaCl, 1mM EDTA, 0.05M Glycine.
Protein was determined by OD₂₈₀ absorbance.

Concentration: 1mg/ml

Physical State: Frozen liquid

Purity: > 95% based on SDS-PAGE

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS

Research Applications

ELISA: to be tested at 0.1~1ug/ml

Coating DM1-antibody	1ug/well				
Purified antibody dilution(ug/ml)	1	0.5	0.2	0.1	0
1A6(1mg/ml)	>3.5	>3.5	>3.5	>3.5	0.047
