



京天成生物技术（北京）有限公司
AbMax Biotechnology Co., LTD
北京经济技术开发区科创十四街 99 号 18 号楼
2 单元 201
99 Kechuang 14th St, Bldg 18, 2-201
BDA, Beijing, China 101111
Tel: 86-10-59755729, Fax: 86-10-59755718
E-mail: info@antibodychina.com
www.abmaxcn.com

Mouse anti-Adalimumab(Humira) 11D5,
Antigen non-competing
(mouse monoclonal IgG)
Catalog # 05-0099

Immunogen: Adalimumab(Humira).

Specificity: Recognizes Adalimumab(Humira) specifically,
no cross reactivity with other humanized antibodies.

Formulation: 0.1mg of mouse monoclonal IgG in 0.1ml
of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02%NaN₃.

Protein was determined by OD₂₈₀ absorbance.


Physical State: Frozen liquid

Purity: > 95% based on SDS-PAGE

Storage: Store the product at -20 °C .
Product is stable for about 6 weeks at
2-8 °C as an undiluted liquid. Prepare
working dilution prior to use. Avoid
repeated freezing and thawing.

FOR IN VITRO RESEARCH USE ONLY

Research Applications



10—5000ng/ml;
LOD: 10ng/ml.