



京天成生物技术（北京）有限公司  
AbMax Biotechnology Co., LTD  
北京经济技术开发区科创十四街 99 号 18 号楼  
2 单元 201  
99 Kechuang 14th St, Bldg 18, 2-201  
BDA, Beijing, China 101111  
Tel: 86-10-59755729, Fax: 86-10-59755718  
E-mail: info@antibodychina.com  
www.antibodychina.com

## Mouse anti-Panitumumab (Vectibix) 2B10

(mouse monoclonal IgG)

Catalog # 05-0121B

**Immunogen:** Panitumumab (Vectibix)

**Specificity:** Recognizes Panitumumab (Vectibix) specifically,  
no cross reactivity with other humanized antibodies.

**Formulation:** 1mg of mouse monoclonal IgG in 1ml  
of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02%NaN<sub>3</sub>.

Protein was determined by OD<sub>280</sub> absorbance.

**Physical State:** Frozen liquid

**Purity:** > 95% based on SDS-PAGE

**Storage:** Store the product at -20 °C .  
Product is stable for about 6 weeks at  
2-8 °C as an undiluted liquid. Prepare  
working dilution prior to use. Avoid  
repeated freezing and thawing.

### FOR IN VITRO RESEARCH USE ONLY

#### Research Applications

##### Indirect ELISA:

	Purified antibody (ug/ml)					
Dilution	1	0.5	0.05	0.005	0.0005	PC/NC
1st antibody						
2B10	999	999	999	999	1.227	0.37