



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-Rabies Virus, Clone 1C7

(mouse monoclonal IgG), FITC-conjugated

Catalog # 12-0020

Antibody Class: Mouse monoclonal IgG produced *in vitro* by mouse-mouse hybridoma 1C7, purified by Protein G-Sepharose chromatography and cross-linked to FITC.

Immunogen: Rabies virus

Formulation: 100 reactions IgG-FITC conjugate in 0.2ml PBS, 25% glycerol.

Physical State: Frozen liquid

Storage: Store the product away from light at -20°C. Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS

Research Applications

Immunohistochemistry: use at 1:25 dilution.
