



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-Enterovirus 71, Clone 13D7

(mouse monoclonal IgG), HRP-conjugated

Catalog # 12-0003

Lot # 100271

Antibody Class: Mouse monoclonal IgG produced *in vitro* by mouse-mouse hybridoma 13D7, purified by Protein G-Sepharose chromatography and cross-linked to horseradish peroxidase.

Immunogen: Synthetic peptide corresponding to a.a 163-177 of Capsid Protein VP1 of Enterovirus 71.

Formulation: 0.2mg of IgG-HRP conjugate in 0.2ml of PBS, 25% glycerol.

Physical State: Frozen liquid

Storage: Store the product away from light at -20°C. Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

**FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS**

Research Applications

ELISA: use at 1:500 to 1:1000
