



京天成生物技术（北京）有限公司  
AbMax Biopharmaceuticals Co., LTD  
北京经济技术开发区科创十四街 99 号 18 号楼  
2 单元 201  
99 Kechuang 14th St, Bldg 18, 2-201  
BDA, Beijing, China 101111  
Tel: 86-10-59755729, Fax: 86-10-59755718  
E-mail: info@antibodychina.com  
www.antibodychina.com

**MOUSE ANTI-PNEUMOCOCCAL  
POLYSACCHARIDE(9N)  
(mouse monoclonal IgG)**

**Catalog # 05-0119B**

**Immunogen:** The 9N Serotype of Pneumococcal capsular polysaccharide (PCPS).

**Specificity:** Recognizes the 9N serotype of PCPS specifically, no cross reactivity with other serotypes of PCPS (1、3、5、4、6A、6B、7F、14、18、19A、19F、23F,etc).

**Formulation:** 1mg of mouse monoclonal IgG in 1ml of 0.1M Tris-HCl, 0.15M NaCl, 0.1mM EDTA, 0.05M Glycine. Protein was determined by OD<sub>280</sub> absorbance.

**Concentration:** 1mg/ml

**Physical State:** Frozen liquid

**Storage:** Store the product at -20°C. Product is stable for about 6 weeks at 2-8°C as an undiluted liquid. Prepare working dilution prior to use. Avoid repeated freezing and thawing.

**Purity:** > 95% based on SDS-PAGE

**FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS**

---

**Research Applications**

**Indirect ELISA:** to be used at 0.5~2µg/ml.

**Immunoturbidimetry:** to be used at 1:2~ 3 dilutions.

---